

Ultrasonic Cleaners





LU1



LU2 & LU3



LU4



LU5



LU6



LU7

Specifications and Ordering Information

Model No.	LU1	LU2	LU3	LU4	LU5	LU6	LU7
Material	Stainless steel including cover				Stainless steel including cover and drain		
Power Supply	110 or 220V, 50/60 Hz						
Power	50 Watt Ultrasonic + 100 Watt Heating	80 Watt Ultrasonic + 100 Watt Heating	50 Watt Ultrasonic + Heating	100 Watt Ultrasonic + 200 Watt Heating	200 Watt Ultrasonic + 100 Watt Heating	500 Watt Ultrasonic + Heating Function	300 Watt Ultrasonic + Heating function
Number of Transducer/ Transducer Frequency (KHz)	1/40	1/28	1 x 40	2 x 40	4 x 40	10 x 40	6 x 40
Tank Size (L)	0.7	2	2	4.5	6	20	10
Tank Dimension (L x D x H) mm	176 x 109 x 65	176 x 164 x 100		327 x 176 x 100	327 x 176 x 150	353 x 327 x 150	327 x 265 x 150
Unit Dimension (L x D x H) mm	176 x 111 x 170	176 x 185 x 212		176 x 190 x 240	176 x 190 x 280	430 x 340 x 430	390 x 290 x 320
Unit Weight (kg)	1.9	2.3		4.5	5.5	13.5	8.5
Timer (min)	0 - 30						
Heating Temp. (°C)	Max. 70						
Catalog No.	37950101	37950102	37950103	37950104	37950105	37950106	37950107

Accessories

Item	Catalog No.
Stainless steel basket	37950111
Rack (Ring holder)	37950112